

Approach for Determining the Background Library Spectra Using a Pulsed D-T Source

Vincent DiNova
NC State University

2018 CNEC Workshop
February 8-9, 2018
Raleigh, North Carolina

Table of Contents

- Monte Carlo Library Least Squares (MCLLS)
- Previous methods of determining background
- Pulsed Neutron Generator (PNG) die-off time
- Determining background using PNG
- KSU/RDRS equipment capabilities
- Next Steps
- Discussion/Questions



Monte Carlo Library Least Squares

- Same principle as least squares except library spectra are used as the input
- $y_i = \sum_{j=1}^m x_j a_{ij} + E_i, i = 1, n$
 - y_i is the counts per channel i
 - x_j are linear coefficients for each element j
 - a_{ij} are the library spectra, or counts in channel i of element j
 - E_i is random error in counts in channel i
- One determination missing is the use of a background library spectra

Previous Methods for Determining Background

- Take background counts with a separate tool, or on a separate run
 - Neither option cost effective
- Residual analysis
 - Perform original MCLLS fitting
 - Using residuals, determine missing elements/libraries
 - Add elements/libraries to MCLLS for final fitting
 - Not a true on-line analysis technique
- Generate “best guess” background composition libraries
 - Using individual nuclides
 - Adds complexity to MCLLS calculation
 - Using semi-empirical background fitting functions

Pulsed Neutron Generator (PNG) die-off time



- Pulse from PNG
- Increase in gamma rays from inelastic scatter and neutron capture
- After 15-20 microseconds, inelastic gamma rays no longer contribute to total spectrum
- Neutron capture gamma rays continue exponential decrease until next pulse

Determining Background Using PNG

- Utilizing complex timing equipment, multiple spectra can be generated during the same pulsing sequence
- The burst spectrum is taken while the PNG is emitting fast neutrons, and includes both capture and inelastic gamma rays
- The capture spectrum is taken just after the PNG is deactivated, when only neutron capture gammas are being emitted
- After repetitive pulsing over a few milliseconds, a background spectrum can be collected before beginning the pulsing sequence again

KSU/RDRS equipment capabilities



- CAEN desktop digitizer
 - Model DT5730B
- Thermofisher B320 D-T Generator
- Canberra 2007P NaI(Tl) detectors
 - 3"x3"
- LND 252 He-3 detectors
 - 1"x8"



Next Steps

- Run preliminary tests to generate data
 - Very large data files
 - Determine best method to parse data into usable form that is also sharable
- Optimize timing sequence to gather only pertinent information in each timing segment
- Apply technique to full D-T cases at KSU testing facility

Questions/Discussions?

- CAEN. *User Manual UM3148*. (2016).
- Gardner, Robin P., Xu, Libai, “Status of the Monte Carlo Library Least-Squares (MCLLS) Approach for Non-Linear Radiation Analyzer Problems.” *Radiation Physics and Chemistry*. 78. (2009).
- Mildenberger, Frank, Mauerhofer, Eric, “Cyclic Neutron Activation Analysis of Large Samples with a Pulsed 14 MeV Neutron Source.” *Journal of Radioanalytical and Nuclear Chemistry*. 311. (2017).
- Tan, H. et al., “A Digital Spectrometer Approach to Obtaining Multiple Time-Resolved Gamma-Ray Spectra for Pulsed Spectroscopy.” *Nuclear Instruments and Methods in Physics Research B*. 263. (2007).
- Tan, H. et al., “A Multiple Time-Gated System for Pulsed Digital Gamma-Ray Spectroscopy.” *Journal of Radioanalytical and Nuclear Chemistry*. Vol. 276, No. 3. (2008).



This work was sponsored in part by the NNSA
Office of Defense Nuclear Nonproliferation R&D
through the Consortium for Nonproliferation Enabling Capabilities

