


RDRS Thrust Area Update

Walter McNeil


Mechanical & Nuclear Engineering Dept.
Kansas State University
Manhattan, KS 66506

June 13, 2018



- The DT generator beam current has been unlocked with help from Fabian
 - DT generator can now be operated from 20-60 μA beam current,
 - Voltage range is 40-80 kV
- The best combination of current and voltage is 40 μA :40kV
 - Lowest dead time and noise with stable neutron output
- Generator noise is now entirely removed from spectrum using decimation settings

We only assume the time response of the neutron output follows:




- Neutron detectors have significantly reduced dead time thanks to reduction of output

Properties

[Spectrum](#) [Board Info](#)

Board Model: DT5730B
 Serial Number: 1492
 Channel: 0
 ADC Channels: 1024
 Real Time: 134.929 s
 Live Time: 133.614 s
 Dead Time: 0.974 %
 Spect. Counts: 141692
 Total Counts: 71174
 ICR: 5.275E+002 Cnt/s
 OCR: 1.050E+003 Cnt/s
 Histogram: BF3_Near




- Neutron detectors have significantly reduced dead time thanks to reduction of output

Properties

[Spectrum](#) [Board Info](#)

Board Model: DT5730B
 Serial Number: 1492
 Channel: 2
 ADC Channels: 1024
 Real Time: 108.931 s
 Live Time: 108.495 s
 Dead Time: 0.400 %
 Spect. Counts: 12878
 Total Counts: 3232
 ICR: 2.967E+001 Cnt/s
 OCR: 1.182E+002 Cnt/s
 Histogram: 3He_Far



- NaI detectors also behaving well, similarly low dead time

Properties ⌵ ✕


[Spectrum](#) [Board Info](#)

Board Model: DT5730B
 Serial Number: 1492
 Channel: 1
 ADC Channels: 2048
 Real Time: 207.930 s
 Live Time: 203.915 s
 Dead Time: 1.931 %
 Spect. Counts: 808845
 Total Counts: 740128
 ICR: 3.560E+003 Cnt/s
 OCR: 3.890E+003 Cnt/s
 Histogram: NaI_Near

Properties




[Spectrum](#) [Board Info](#)

Board Model: DT5730B
 Serial Number: 1492
 Channel: 3
 ADC Channels: 2048
 Real Time: 179.931 s
 Live Time: 178.609 s
 Dead Time: 0.735 %
 Spect. Counts: 89548
 Total Counts: 88114
 ICR: 4.897E+002 Cnt/s
 OCR: 4.977E+002 Cnt/s
 Histogram: NaI_Far



External Timing Trigger

- DT Generator is now connected to a function generator
- Switching noise is synchronized with square wave



- List-Mode data provides timestamps ($2\mu\text{s}$ intervals) and channel data (15 bit) for each pulse event in three numbers
- ASCII format is easy to interpret, but file sizes are large

TIMESTAMP CHANNEL EVENT_TAG

Sensor on Ch0

```
5_1_2018-Test_ch000 - Notepad
File Edit Format View Help
HEADER0:1281
HEADER1:1792
HEADER2:513
HEADER3:770
HEADER4:35588
793462 6072 0
2487657 4630 0
2723981 6393 0
2973176 5480 0
3576909 3221 0
3806627 6007 0
4015208 4546 0
4056108 -32768 0
4057038 -32768 0
4264107 2689 0
4564343 5957 0
4765540 5563 0
4860340 4674 0
5036171 5093 0
5747110 5924 0
5785175 3970 0
7409749 5579 0
9065690 5932 0
9140009 5665 0
9217795 5768 0
9588844 5500 0
9700273 6798 0
10001978 5710 0
10428728 4630 0
```

Pulse Driver on Ch4

```
5_1_2018-Test_ch004 - Notepad
File Edit Format View Help
HEADER0:1281
HEADER1:1792
HEADER2:513
HEADER3:770
HEADER4:35588
374113 1991 0
874097 1991 0
1374082 1991 0
1874066 1991 0
2374051 1990 0
2874035 1991 0
3374020 1991 0
3874004 1991 0
4373989 1989 0
4873973 1989 0
5373958 1990 0
5873942 1991 0
6373927 1990 0
6873912 1990 0
7373896 1990 0
7873881 1991 0
8373865 1991 0
8873850 1990 0
9373834 1991 0
9873819 1990 0
10373803 1990 0
10873788 1990 0
11373772 1991 0
11873757 1990 0
```


- The DT generator beam current has been unlocked with help from Fabian
 - DT generator can now be operated from 20-60 μA beam current,
 - Voltage range is 40-80 kV
- The “**best**” combination of current and voltage is 40 μA :40kV
 - Low dead time and noise with **good neutron** output
- Noise is now entirely removed from spectrum using decimation settings

I question this now?





Three Tap Water Experiments

06/11/18



Long Vo

Alan Reinke & Alexi Rojas

Repeatability of time-integrated spectra. NaI



Repeatability of Time-Integrated Spectra (BF3 and He-3)



Dead-time is acceptable

- TapWater00

1. NearNaI: 1.2%
2. FarNaI: 0.6%
3. NearBF3: 9.9%
4. FarHe3: 12%

- TapWater01

1. NearNaI: 1.2%
2. FarNaI: 0.6%
3. NearBF3: 13%
4. FarHe3: 9.6%

- TapWater03


1. NearNaI: 1.2%
2. FarNaI: 0.6%
3. NearBF3: 12%
4. FarHe3: 13%

$$\frac{RealTime - LiveTime}{RealTime} \times 100$$

Time-Dependent Spectral Data


- Time progression of pulse energy spectra

Calibration NAI near



From TapWater00


- Time progression of energy spectra



From TapWater00



Red dots identify error flags on time-stamp
(likely pile-up, or unusual waveform)

- Time progression of energy spectra




From TapWater00

- Time progression of energy spectra over 2 seconds




Time-Dependent Gross Count-Rate


TW00: NearBF3:1st4sec




TW00: FarNal: 1st4sec



TW00: NearNal: 1st4sec




TW00: FarHe3: 1st4sec




Cyclical Output Behavior


1st 4 sec
1646cts




2nd 4sec
1737cts



about 15 minutes in
1730cts







about 35 minutes in
1796cts




A Look at Die-Away

- Time between every trigger is 1000 micro-second
- Binning Time after every trigger into 10us bins
- Y-axis is count per second
- A clock rollover is ~4295triggers or $4295 \times (1000\text{us}) = 4.295\text{second}$


Die Away




- What energies are we losing because of cut off?
- Looking at regions (20us-140us), (140us-500us), and (500us-1000us) in NaI_Near



From TapWater00



From TapWater01



Fast Neutron and Gamma Burst measured with
EJ-299-34 Plastic and SiPM on THERMOFISHER
A325 DT Generator 2KHz and 20% duty cycle




Figure 6: Fast neutron and gamma count rates over the course of a 2 kHz, 20% duty cycle burst period.

R. M. Preston, "Neutron generator burst timing measured using a pulse shape discrimination plastic scintillator with silicon photomultiplier readout"